


Supercritical Fluid-based Processing Technologies from IIT Bombay

Supercritical Fluid-based (SCF) Processing has been a significant research focus at the Department of Chemical Engineering at IIT Bombay (IITB) for over three decades. Currently IITB houses a wide range of experimental facilities: from bench-scale setups to pilot plant. The principal areas of research are:

- Extraction of natural products
- Micronization and Encapsulation of high-value nutraceuticals and pharmaceuticals
- Preservation of liquid and solid foods
- Reactions in supercritical fluids
- Pressurized Water-based extraction processes for natural medicinals and nutraceuticals

Key SCF systems available at IITB:

1. **Supercritical Fluid Extraction System:** SCFE is a two-step process (see schematic below) which uses a dense gas as a solvent e.g., carbon dioxide (CO_2) for extraction, above its critical temperature (31°C) and critical pressure (74 bar). The feed, generally ground solid, is charged into the extractor. Supercritical CO_2 is fed to the extractor through a high pressure pump (100-500 bar). The extract laden CO_2 is sent to a separator (60-120 bar) via a pressure reduction valve. At reduced temperature and pressure conditions, the extract precipitates out in the separator. The extract free CO_2 stream, leaving the separator is then recycled to the extractor. In the case of liquid feed, the extractor is modified into a column through which feed and the supercritical CO_2 is fed either co-currently or concurrently.







Supercritical Fluid Extraction Pilot Plant at IITB

2. Supercritical Fluid Micronization/Encapsulation System: Particle design (micronization and encapsulation) is presently a major development of supercritical fluids applications, mainly in the pharmaceutical, nutraceutical, cosmetic and specialty chemistry industries. Typically, the supercritical fluid is used as an anti-solvent that causes precipitation of the substrate(s) dissolved initially in a liquid solvent. The supercritical fluid is pumped to the top of the high pressure vessel (see schematic below) by a high pressure pump while the active substance solution is introduced into the high pressure vessel through a nozzle. To produce small liquid droplets in the nozzle, the liquid solution is pumped at a pressure higher than the vessel operating pressure. Particles precipitate due to the anti-solvent action of the supercritical carbon dioxide, and are collected on a filter at the bottom of the vessel. The fluid mixture (supercritical fluid plus solvent) exits the vessel and flows

depressurization tank where the conditions allow gas-liquid separation.



3. Pressurized Hot Water Extraction System: At normal conditions the dielectric constant value of water is quite high, but at sub- or supercritical conditions, it is lowered water can mimic the behaviour of both polar and apolar solvents, and hence allowing the solubilization of a wide variety of moieties normally not very soluble in water. Subcritical water extraction or Pressurised Hot Water Extraction (PHWE) is based on this unique solvent property of water. A very wide range of herbal products may be extracted efficiently using this process. The PHWE system schematic is shown below.



- V1 = Opening valve of gas reservoir
- V2 = Inlet valve to extractor
- V3 = Sample collection valve
- V4 = Gas Purging valve
- V5 = needle valve for collection of extracted soln

Hot water is pressurized using gaseous CO₂ at a moderate (2–20 bar) pressure. The pressurized water is then used to extract the analyte from the plant matrix. The extraction

efficiency depends upon: processing temperature and pressure, particle size of plant matrix, water to plant matrix ratio, etc. PHWE is suitable for a extraction of a various natural molecules: flavonoids, essential oils, saponins, etc.

Services offered to the industry:

The IIT, Bombay research team undertakes both fundamental and applied research in diverse areas of supercritical fluid technology. Apart from the study of a large variety of specific systems requiring SCF processing, the research experience has been translated into development of *process* and *engineering designs* of SCF systems (bench top to commercial scale) for:

- (i) *SCF-based extraction of natural products (spices, medicinal herbs, natural colours, flavours and fragrances, etc.)*
- (ii) *SCF-based micronization and encapsulation of nutraceuticals and pharmaceuticals*
- (iii) *Pressurized hot water extraction of herbals*

The specific services offered to the industry are:

- *Technology Consultation: Process Optimization; Selection / Choice of Viable Products; Selection of Optimum Process Plant Configuration, Feasibility Analysis*
- *Generation of Extracts for Test Marketing*
- *Contract Research*

Contact:

The Principal R&D Team Members from IIT Bombay are faculty with the Department of Chemical Engineering

- Prof Sandip Roy (sr@che.iitb.ac.in, 91-22-2576 7249)
- Prof Madhu Vinjamur (madhu@che.iitb.ac.in, 91-22-2576 7218)
- Prof Mamata Mukhopadhyay (mm@che.iitb.ac.in; 91-22-2576 7248)
- For more information about Supercritical Fluid Processing Technologies from IITB please visit: <http://www.che.iitb.ac.in/online/services/technology-transfers>
- Email: scfe@che.iitb.ac.in